ENGLISH

INSTRUCTIONS FOR USE– DISPOSABLE CAUTERY

This cautery is disposable and is alkaline-battery powered, supplied **STERILE** and intended for single patient use. **DO NOT** resterilize this device. Contents are sterile if the package has not been opened or is not damaged. For questions or additional information on our complete line of cautery's and electrosurgical devices, please contact Kirwan Surgical Products LLC (KSP) at (781) 834-9500 or toll free at 1-888-547-9267, or on the web at www.ksp.com.

INTENDED USE: Cautery units are intended for general surgical use to cauterize small blood vessels in order to help reduce blood loss and keep the site clear of excess blood.

CAUTION: Any use of this cautery for tasks other than for which it is indicated, will usually result in a damaged or broken cautery.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

CAUTION: As with all alkaline batteries, store this device at non-elevated temperatures to avoid battery power loss.

SETUP and USE

A protective cap covers the switch and tip. Carefully remove cap and activate switch to heat Cautery tip. Replace cap when not in use. At low temperatures, it glows a dull gray at the tip; at high temperatures, it glows a bright red. For variable Cautery version, turn dial to right to increase tip temperature to desired level.

PROPER DISPOSAL of cautery's possibly contaminated with blood, tissue, or other potentially infectious material present a biological risk and must be discarded in a closable, leak-proof, puncture-resistant receptacle, that is adequately labeled (e.g., color coding or symbology) for easy identification as biohazard waste.

WARNING: Do not reuse or re-sterilize. Patient cross contamination and/or product malfunction may occur.

DANGER FIRE RISK WARNINGS


Replace cap when not in use.

Do not use this cautery in the presence of flammable gasses or materials such as facial hair, preparation agents, alcohol vapors, drapes, or gowns, or oxygen rich environments.

BEFORE DISCARDING break off and completely remove the wire tip with hemostat or scissors and replace safety cap.

Step 1 – Break off tip

Step 2 – Replace cap to immobilize Actuator Button

Failure to follow these instructions could result in fire.

REF Catalog Number
LOT Batch Code
 Use by date
 Do not use if package is damaged
 Do not reuse
 Do not resterilize

Consult Instruction for Use
NO LATEX Not made with natural rubber latex.
STERILE R Sterile -Sterilized using irradiation
 Warning, Cautions and Precautions

Rx Only CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.
Product of USA

Kirwan Surgical Products LLC
180 Enterprise Drive
Marshfield, MA 02050, USA
Tel.: (781)-834-9500

Medical Product Service GmbH
Borngasse 20
35619 Braunfels, Germany

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